

LABORHITS

pH-Meter Seven2Go™ S2

- One-hand operation
- Simple menu structure
- Waterproof (IP67)



Items supplied:

- S2-Meter:** Instrument, electrode clip, operating instructions, quick guide, test certificate, declaration of conformity and batteries.
- S2-Standard kit:** As S2-Meter, but also with InLab Expert Pro-ISM-IP67, buffer sachets for pH 4.01, 7.00, 9.21 and 10.00.
- S2-Field kit:** As S2-Standard kit, but also with bottles for calibration and uGo™ carrying case.
- S2-Food kit:** As S2-Field kit, with InLab Solids Pro-IP67 instead of InLab Expert Pro-ISM-IP67.
- S2-Light kit:** As S2-Standard kit, with InLab Versatile Pro instead of InLab Expert Pro-ISM-IP67.

Specifications

pH	Measuring range:	-2 to 20 pH
	Resolution:	0.01/0.1
	Accuracy:	±0.01
mV	Measuring range:	-1999 ... +1999 mV
	Resolution:	1
	Accuracy:	±1
	Relative mV-measurement:	yes
Temperature	Measuring range:	-5.0 ... 105.0 °C
	Resolution:	0.1 °C
	Accuracy:	±0.5 °C

Advantages:

- Portable routine model for pH measurement
- New, intuitive menu structure with fast and easy T-Pad navigation in 12 languages
- Robust design with IP67 protection
- Intelligent Sensor Management (ISM®) functionality reduces the risk of errors
- Automatic buffer recognition and temperature compensation
- 4 predefined buffer groups and 1 user defined buffer group
- Automatic, manual and time dependant endpoint recording
- Data memory storage for 200 data points
- Power requirements: 4 x 1.5V AA or 1.3V NiMH accumulators

Typ	Cat. No.
S2-Meter	6.265 414
S2-Standard kit	6.265 415
S2-Field kit	6.265 416
S2-Food kit	6.265 417
S2-Light kit	6.265 418