


LABORHITS

Multi-Purpose-Rack

For 5 ml, 15 ml and 16 mm Tubes

Height Adjustable

Autoclavable



This multi-purpose rack is designed with the flexibility to hold 20 x 5 ml and 20 x 15 ml tubes or to hold 40 x 5 ml or 40 x 15 ml tubes.

In addition, the rack can also hold clinical tubes of a similar diameter (16 mm) and height. The rack incorporates two moving hinged upper shelves that can be rotated to an up or down position to change height to accommodate the size of the tube being used. Unique tear drop shaped openings (top tray) for better retention of tubes.



Ergonomic handles for ease of use. Imprinted grid on both sides of the upper shelves for quick identification and retrieval.

Dimensions (L x W x H): 309 x 107 x 77 mm

Array	For tubes diam. mm	Colour	PK	Cat. No.
40	16	green/blue	1	6.281 157
40	16	purple/natural	1	6.281 249