


# LABORHITS

## WHEATON® Bioreactors CELLine™



- Disposable and ready-to-use
- High cell density and high product concentration
- No additional equipment required for operation

New line of membrane driven flasks for high density cell culture. The flasks are designed to enhance small scale bio-production for antibodies and proteins. Conventional *in vivo* or *in vitro* cell culture methods can be laborious, may have low cell density and require significant purification. CELLine™ flasks address these three areas of limitation observed in static tissue culture flasks.



Type	Culture type	Media Compartment ml	Cell Compartment ml	PK	Cat. No.
CELLine™ 1000	Suspension	1000	15	1	4.658 779
CELLine™ 1000	Suspension	1000	15	1	4.658 780
CELLine™ 1000-AD	Adherent	1000	15	1	4.658 781
CELLine™ 1000-AD	Adherent	1000	15	1	4.658 782
CELLine™ 350	Suspension	350	5	1	4.658 777
CELLine™ 350	Suspension	350	5	1	4.658 778