

LABORHITS **PROMO**

Techne ³PRIMEBASE Thermal Cycler

**Promotion until
30.06.2016**

Fast heating rate 3 °C/sec

Easy intuitive operation

Uniform temperature distribution

The ³PRIMEBASE Thermal Cyclers provide easy programming and feature a full colour touchscreen. Pre-installed templates allow the user to quickly create even the most complicated programmes via the intuitive touch interface. All Prime thermal cyclers have the ability to transfer programmes between instruments via USB or to your PC. The internal memory can hold 1000 programmes with unlimited storage available on USB data stick.

Specifications

Sample capacity: 24 x 0.2 ml, 18 x 0.5 ml

Maximum heating rate: 3 °C/sec

Maximum number of programs stored: 1000

Block temperature range: 10 °C – 100 °C

Block uniformity at 55 °C: ±0.3 °C

Block options: 24 x 0.2 ml

(compatible with 8 tube strips) and 18 x 0.5 ml

Temperature accuracy at 55 °C: ±0.25 °C

Adjustable heated lid temperature: 35 °C to 115 °C or off

Program interface: 3.5 inch touchscreen

Dimensions (L x W x H):

350 x 210 x 180 mm

Weight: 6 kg

Voltage: 100-230 V, 50-60 Hz

Power: 155 W



**Save
€ 600**



Key Features:

- Fast ramp rate up to 3.0 °C/s
- High temperature accuracy (±0.25 °C and at 55 °C) and block uniformity (±0.3 °C and at 55 °C)
- Colour touchscreen user interface for fast program setup
- Data transfer via USB
- Techne 4 year warranty

Description	Cat. No.	List price	Offer price
³ PRIMEBASE Thermal Cycler, 24 x 0.2 ml	9.595 502	2,500-€	1,895 €
³ PRIMEBASE Thermal Cycler, 18 x 0.5 ml	9.595 501	2,500-€	1,895 €

Valid until 30.06.2016 and excl. of VAT.