

LABORHITS

Micro Dispenser VITLAB[®] piccolo 1



Serial pipetting of the innovative kind

For dispensing tiny quantities

Rotates freely 360°

- connecting thread GL 28
- precise dosing
- one hand operation
- low weight
- dispensing tube can be rotated 360°
- only high-quality materials, such as PTFE, PFA, ETFE, FEP, borosilicate glass, and platinum-iridium come in contact with the media

For dispensing small quantities of liquids in all areas of biochemical and medical research.

Even the smallest quantities can be dispensed directly from the bottle with the VITLAB[®] piccolo 1 - helpful particularly for serial dispensing operations.

The ergonomic design makes dispensing effortless and stress-free. The VITLAB[®] piccolo 1 can be operated with only one hand. Use the thumb depress the volume dispensing button, just as with a pipette,

and a reset mechanism refills the volume automatically.

The dispensing tube can be rotated over 360° so that it is always optimally situated with respect to the bottle label.

Included in delivery: VITLAB[®] piccolo 1 with GL 28 connecting threads and mounting tool. Also available with DAkKS calibration certificate or individual quality certificate (on request).

Capacity µl	Accuracy max. vol. (≤ ± R%)	Precision max. vol. (≤ CV%)	PK	Cat. no.
100	3,0	0,4	1	9.284 045
200	2,5	0,4	1	9.284 044
250	2,0	0,4	1	9.284 046
500	1,5	0,3	1	9.284 047
1000	1,0	0,2	1	9.284 048