

LABORHITS

Cryogenic Boxes true north[®] Flatpack, PP

Safe storage up to -80 °C



Easy self-assemble boxes shipped flat-packed saving space in the laboratory.

- Made of corrugated polypropylene material
- Not susceptible to moisture
- Reusable
- Dividers included with the boxes

For tubes ml	Colour	Dimensions (W x D x H) mm	Array	PK	Cat. No.
0.2	white	133 x 130 x 30	144	10	6.281 163
0.2	blue	133 x 130 x 30	144	10	6.281 164
0.2	red	133 x 130 x 30	144	10	6.281 165
0.5	white	133 x 130 x 42	81	10	6.281 166
0.5	blue	133 x 130 x 42	81	10	6.281 167
0.5	red	133 x 130 x 42	81	10	6.281 168
1.5/2.0	white	133 x 130 x 51	81	10	6.281 158
1.5/2.0	blue	133 x 130 x 51	81	10	6.281 161
1.5/2.0	red	133 x 130 x 51	81	10	6.281 162
5.0	white	133 x 130 x 75	25	10	6.281 169
5.0	blue	133 x 130 x 75	25	10	6.281 170
5.0	red	133 x 130 x 75	25	10	6.281 171
15.0	white	145 x 147 x 120	36	10	6.281 172
15.0	blue	145 x 147 x 120	36	10	6.281 173
15.0	red	145 x 147 x 120	36	10	6.281 174
50.0	white	145 x 147 x 120	16	10	6.281 159
50.0	blue	145 x 147 x 120	16	10	6.281 160
50.0	red	145 x 147 x 120	16	10	6.281 536