

LABORHITS

Wireless module IDS WLM-S/WLM Kit

Same module for all parameters

Wireless measurement

Securely connected

**Complete set:
Wireless module IDS WLM KIT:**



No cables, no tugging, full flexibility in all areas where cables or meters would disturb, even under laboratory hoods or laminar flow benches.

Suitable for IDS portable meters:
MultiLine® Multi 3630 IDS / 3620 IDS / 3510 IDS
Suitable for IDS laboratory meters:
inoLab® Multi 9630 IDS / 9620 IDS / 9310 IDS



Specifications

Range: About 10 m

Continuous operating times

IDS pH/ORP electrodes: 60 h

IDS Conductivity cells: 30 h (conductivity up to 1 S)

IDS Oxygen sensor: 9 h

IDS Turbidity sensor: 5 h

Wireless transmission: Bluetooth® LE

Signal: LED RGB

Protection class: IP43

Power supply: Re-chargeable LiPo battery 230 mAh

Type	Description	PK	Cat. No.
Wireless module IDS WLM-S	for plug head sensor	1	6.283 617
Wireless module IDS WLM KIT	incl. USB-Charger, USB-Power supply	1	6.283 612
WLM-Charger	USB-Charger, cascable	1	6.283 618