

# LABORHITS

## VACUU-VIEW and VACUU-VIEW extended

### Vacuum measuring instruments

Rough vacuum gauge  
**VACUU-VIEW** -  
Measuring with  
highest precision from  
atmosphere down to  
0.1 mbar

Compact design with  
integrated sensors

Compatible with gauge  
DCP 3000 and controller  
CVC 3000

Rough and fine vacuum gauge  
**VACUU-VIEW extended**  
with an extended measurement  
range from atmosphere down to  
 $10^{-3}$  mbar

Easy to read,  
illuminated display

Easy, menu driven  
handling

Chemically resistant  
and robust



Specifications	VACUU-VIEW	VACUU-VIEW extended
Measuring range	1100 - 0.1 mbar (hPa)	1100 - 0.001 mbar (hPa)
Measuring range	825 - 0.075 torr	825 - 0.001 torr
Measuring principle	ceramic diaphragm sensor with capacitive readout	ceramic diaphragm sensor + ceramic jacketed Pirani sensor
Accuracy of measurement	±1 mbar (hPa)	±15 % of indicated value in the range from 10 to 0.01 mbar (hPa) ±3 mbar for >10 mbar (hPa)
Vacuum connection	KF DN 16 / hose nozzle DN 6/10 mm	
Dimensions (L x W x H)	103 x 62 x 50 mm	
Weight	190 g	
Power supply	100 - 230 V, 50/60 Hz (CEE/CH/UK/AUS/CN)	
Cat. No.	6.268 865	6.268 864

